


石英NDフィルター 5%



波長 (nm)	RawData ...
300.00	5.766
320.00	6.470
340.00	6.849
360.00	6.931
380.00	6.848
400.00	6.660
420.00	6.455
440.00	6.242
460.00	6.063
480.00	5.926
500.00	5.846
520.00	5.801
540.00	5.767
560.00	5.743
580.00	5.708
600.00	5.677
620.00	5.644
640.00	5.613
660.00	5.566
680.00	5.542
700.00	5.515
720.00	5.506
740.00	5.486
760.00	5.464
780.00	5.513
800.00	5.529
820.00	5.551
840.00	5.580
860.00	5.604
880.00	5.668
900.00	5.717
920.00	5.723
940.00	5.750
960.00	5.779
980.00	5.806
1000.00	5.825
1020.00	5.842
1040.00	5.855
1060.00	5.866
1080.00	5.875
1100.00	5.882
1120.00	5.882
1140.00	5.882
1160.00	5.881
1180.00	5.875
1200.00	5.871