



[測定プロパティ]  
 波長範囲(nm): 400.00 ~ 700.00  
 スキャン スピード: 中速  
 サンプルングピッチ: 0.2  
 オート サンプルングピッチ: 有効  
 測定モード: オート

[装置プロパティ]  
 装置の種類: UV-3600Plus シリーズ  
 測光値: 透過率  
 スリット幅: (20)  
 時定数: 0.1 秒  
 光源: 自動  
 光源切替波長: 282.00 nm  
 検出器ユニット: 外部(3検出器)  
 検出器切替波長: 780.00 nm 1800.00 nm

波長 (nm)	0 ...	50 ...	100 ...
400.00	91.7669	64.1286	44.4632
405.00	91.8028	64.2344	44.5612
410.00	91.7640	64.3234	44.6801
415.00	91.7406	64.4485	44.7831
420.00	91.7817	64.5341	44.8913
425.00	91.7796	64.6131	44.9946
430.00	91.7818	64.7226	45.0863
435.00	91.7748	64.7961	45.1714
440.00	91.7443	64.8928	45.2826
445.00	91.7730	64.9608	45.3580
450.00	91.7798	65.0722	45.4560
455.00	91.7952	65.1305	45.5548
460.00	91.8110	65.1869	45.6603
465.00	91.7959	65.2808	45.7252
470.00	91.7628	65.3622	45.7986
475.00	91.7925	65.4165	45.8746
480.00	91.8080	65.5395	45.9696
485.00	91.8164	65.5972	46.0720
490.00	91.8236	65.6912	46.1530
495.00	91.8279	65.7422	46.2297
500.00	91.8531	65.8192	46.3019
505.00	91.8232	65.9118	46.3833
510.00	91.8839	65.9880	46.4743
515.00	91.8877	66.0368	46.5344
520.00	91.8656	66.1156	46.6159
525.00	91.8257	66.1724	46.7090
530.00	91.8643	66.2420	46.8086
535.00	91.8821	66.3297	46.8631
540.00	91.8698	66.3781	46.9272
545.00	91.8836	66.4097	46.9863
550.00	91.8678	66.4962	47.0407
555.00	91.9055	66.5616	47.1416
560.00	91.9051	66.6175	47.2012
565.00	91.9244	66.6571	47.2674
570.00	91.9052	66.7171	47.2937
575.00	91.9128	66.7633	47.3772
580.00	91.8588	66.8462	47.4341
585.00	91.8613	66.9086	47.4798
590.00	91.8739	66.9457	47.5435
595.00	91.8474	66.9848	47.5925
600.00	91.8663	67.0385	47.6454
605.00	91.8717	67.0485	47.7117
610.00	91.8720	67.1298	47.7759
615.00	91.8811	67.1458	47.8290
620.00	91.8853	67.1975	47.8921
625.00	91.8729	67.2529	47.9215
630.00	91.8673	67.3096	47.9955
635.00	91.9024	67.3780	48.0791
640.00	91.8522	67.4274	48.1528
645.00	91.8750	67.4775	48.1859
650.00	91.8928	67.5120	48.2623
655.00	91.9117	67.5854	48.3603
660.00	91.8638	67.6037	48.3409
665.00	91.8351	67.6325	48.3651
670.00	91.8320	67.7040	48.3914
675.00	91.8489	67.7299	48.4950
680.00	91.8204	67.7722	48.5771
685.00	91.8967	67.8560	48.6151
690.00	91.8756	67.8773	48.6904
695.00	91.9000	67.9672	48.7133
700.00	91.7274	67.9990	48.8066