


NDフィルタ- 40%



波長 (nm)	RawData ...
300.00	39.695
320.00	40.086
340.00	40.357
360.00	40.534
380.00	40.501
400.00	40.352
420.00	40.178
440.00	40.006
460.00	39.825
480.00	39.739
500.00	39.766
520.00	39.897
540.00	40.076
560.00	40.281
580.00	40.533
600.00	40.771
620.00	41.019
640.00	41.318
660.00	41.531
680.00	41.797
700.00	42.040
720.00	42.258
740.00	42.524
760.00	42.749
780.00	43.019
800.00	43.261
820.00	43.501
840.00	43.774
860.00	43.994
880.00	44.249
900.00	44.491
920.00	44.717
940.00	44.927
960.00	45.141
980.00	45.355
1000.00	45.562
1020.00	45.736
1040.00	45.907
1060.00	46.087
1080.00	46.243
1100.00	46.400
1120.00	46.540
1140.00	46.679
1160.00	46.816
1180.00	46.944
1200.00	47.063