



波長 (nm)	RawData ...
300.00	0.0063
310.00	0.0039
320.00	0.0023
330.00	0.0041
340.00	0.1480
350.00	7.5669
360.00	34.7929
370.00	62.7570
380.00	79.2881
390.00	87.8143
400.00	91.5227
410.00	92.7404
420.00	93.0775
430.00	93.3332
440.00	93.6658
450.00	94.1476
460.00	94.7102
470.00	95.1783
480.00	95.5257
490.00	95.8479
500.00	96.1103
510.00	96.2226
520.00	96.2952
530.00	96.4505
540.00	96.4318
550.00	96.4754
560.00	96.6326
570.00	96.7674
580.00	96.8852
590.00	96.8469
600.00	96.9752
610.00	97.1427
620.00	97.1974
630.00	97.1168
640.00	96.9637
650.00	96.7845
660.00	96.5567
670.00	96.1765
680.00	95.5928
690.00	94.9495
700.00	94.0269
710.00	93.2736
720.00	92.4188
730.00	91.4643
740.00	90.2000
750.00	89.0147
760.00	87.7394
770.00	86.3828
780.00	84.4878
790.00	82.7244
800.00	81.3782
810.00	79.8619
820.00	78.4843
830.00	77.0021
840.00	75.5824
850.00	74.1510
860.00	72.7657
870.00	71.4199
880.00	70.1043
890.00	68.8089
900.00	67.5474
910.00	66.3719
920.00	65.1360
930.00	64.0042
940.00	62.9640
950.00	61.8799
960.00	60.8641
970.00	59.8783
980.00	58.9587
990.00	58.0949
1000.00	57.2448
1010.00	56.4516
1020.00	55.6598
1030.00	54.9149
1040.00	54.1963
1050.00	53.4887
1060.00	52.8132
1070.00	52.1913
1080.00	51.5999
1090.00	51.0468
1100.00	50.5052
1110.00	50.0126
1120.00	49.5630
1130.00	49.0991
1140.00	48.6675
1150.00	48.2691
1160.00	47.8480
1170.00	47.4724
1180.00	47.1147
1190.00	46.7846
1200.00	46.4568

[測定プロパティ]

波長範囲(nm): 300.00 ~ 1300.00
 スキャン スピード: 中速
 サンプルングピッチ: 0.5
 オート サンプルングピッチ: 無効
 測定モード: オート

[装置プロパティ]

装置の種類: UV-3600Plus シリーズ
 測光値: 透過率
 スリット幅: (20)
 時定数: 0.1 秒
 光源: 自動
 光源切替波長: 282.00 nm
 検出器ユニット: 外部(3検出器)
 検出器切替波長: 780.00 nm 1800.00 nm