



## [測定プロパティ]

波長範囲(nm): 300.00 ~ 1300.00  
 スキャン スピード: 中速  
 サンプルングピッチ: 0.5  
 オート サンプルングピッチ: 無効  
 測定モード: オート

## [装置プロパティ]

装置の種類: UV-3600Plus シリーズ  
 測光値: 透過率  
 スリット幅: (20)  
 時定数: 0.1 秒  
 光源: 自動  
 光源切替波長: 282.00 nm  
 検出器ユニット: 外部(3検出器)  
 検出器切替波長: 780.00 nm 1800.00 nm

波長 (nm)	RawData ...
300.00	0.0064
310.00	0.0240
320.00	0.0253
330.00	0.0253
340.00	0.0206
350.00	0.0822
360.00	4.3553
370.00	23.9147
380.00	47.0870
390.00	64.1119
400.00	74.0654
410.00	79.3873
420.00	81.7441
430.00	83.3355
440.00	84.3723
450.00	85.4278
460.00	86.6687
470.00	86.7242
480.00	87.7320
490.00	88.4312
500.00	88.7241
510.00	88.9852
520.00	89.0492
530.00	89.2521
540.00	89.2513
550.00	89.4316
560.00	89.6047
570.00	90.1391
580.00	89.7113
590.00	90.0492
600.00	90.3925
610.00	90.5814
620.00	90.6329
630.00	90.5380
640.00	90.2045
650.00	89.7467
660.00	89.1962
670.00	87.5604
680.00	86.7225
690.00	85.2672
700.00	83.3563
710.00	81.3597
720.00	79.2283
730.00	77.3653
740.00	75.1090
750.00	72.2348
760.00	69.4806
770.00	67.6669
780.00	60.9164
790.00	50.9718
800.00	49.1640
810.00	47.4983
820.00	45.5532
830.00	43.8812
840.00	42.5448
850.00	40.8513
860.00	39.3324
870.00	37.9833
880.00	36.5592
890.00	35.3252
900.00	34.1054
910.00	32.9895
920.00	31.8550
930.00	30.9545
940.00	29.9662
950.00	29.0681
960.00	28.2138
970.00	27.4312
980.00	26.7013
990.00	26.0290
1000.00	25.4252
1010.00	24.8444
1020.00	24.3135
1030.00	23.8677
1040.00	23.4253
1050.00	23.0016
1060.00	22.6441
1070.00	22.3345
1080.00	22.0229
1090.00	21.7242
1100.00	21.4480
1110.00	21.1882
1120.00	20.9430
1130.00	20.7334
1140.00	20.5681
1150.00	20.4045
1160.00	20.2289
1170.00	20.0481
1180.00	19.9070
1190.00	19.7803
1200.00	19.6938