


LDJ50



[測定プロパティ]

波長範囲(nm): 300.00 ~ 1300.00
 スキャン スピード: 中速
 サンプルングピッチ: 0.5
 オート サンプルングピッチ: 無効
 測定モード: オート

[装置プロパティ]

装置の種類: UV-3600Plus シリーズ
 測光値: 透過率
 スリット幅: (20)
 時定数: 0.1 秒
 光源: 自動
 光源切替波長: 282.00 nm
 検出器ユニット: 外部(3検出器)
 検出器切替波長: 780.00 nm 1800.00 nm

波長 (nm)	RawData ...
300.00	0.0130
310.00	0.0188
320.00	0.0236
330.00	0.0225
340.00	0.0199
350.00	1.7334
360.00	20.6407
370.00	51.3225
380.00	73.0852
390.00	85.0424
400.00	90.6435
410.00	92.8677
420.00	93.6079
430.00	93.8720
440.00	94.1641
450.00	94.3774
460.00	94.5876
470.00	94.8063
480.00	94.8977
490.00	95.2057
500.00	95.3284
510.00	95.4035
520.00	95.3052
530.00	95.5702
540.00	95.6625
550.00	95.7160
560.00	95.9882
570.00	96.0999
580.00	96.4444
590.00	96.7897
600.00	97.0745
610.00	97.2923
620.00	97.3243
630.00	97.3191
640.00	97.3018
650.00	97.0229
660.00	96.5652
670.00	95.9602
680.00	95.5063
690.00	94.4955
700.00	93.4614
710.00	92.1507
720.00	90.5065
730.00	89.5777
740.00	87.9423
750.00	85.7806
760.00	84.6370
770.00	82.8054
780.00	78.1949
790.00	72.4800
800.00	70.7941
810.00	69.0423
820.00	67.4475
830.00	65.6329
840.00	63.9178
850.00	62.4034
860.00	60.6421
870.00	59.2218
880.00	57.7665
890.00	56.4063
900.00	55.0492
910.00	53.7455
920.00	52.5121
930.00	51.2949
940.00	50.1931
950.00	49.1088
960.00	48.0604
970.00	47.0353
980.00	46.0734
990.00	45.1902
1000.00	44.3224
1010.00	43.5439
1020.00	42.7712
1030.00	41.9944
1040.00	41.2956
1050.00	40.6680
1060.00	40.0388
1070.00	39.4413
1080.00	38.9194
1090.00	38.4026
1100.00	37.9074
1110.00	37.4273
1120.00	37.0114
1130.00	36.5621
1140.00	36.1524
1150.00	35.7900
1160.00	35.3798
1170.00	35.0015
1180.00	34.6834
1190.00	34.3526
1200.00	34.0915